



# Operating Manual

# Prosthetic Scanners





# Vinyl Series

A product from smart optics Sensortechnik GmbH



# Legal Notice

Dear customer,

Thank you for choosing a product from smart optics!

If you require additional information or if you have any questions about the device, please contact your smart optics reseller. Our trained resellers around the world provide first-level support as well as training and service for all applications.

We wish you much joy with your scanner!

Your smart optics team

#### PRIZES AND AWARDS

In 2018, Vinyl received the Red Dot Design Award in the category "Industrial Equipment, Machinery and Automation": www.red-dot.org/de/project/vinyl-24269

#### CHANGES AND SERVICE

We reserve the right to implement product changes as part of continuous improvement and technical progress and to make changes to this documentation. You can find the current version of this documentation on our homepage: www.smartoptics.de.

Please note that the current version of this documentation may contain information which does not apply to older devices.

On request we are also glad to send you a printed paper version of the documentation. Please send us an e-mail to communications@smartoptics.de.

If you are looking for a reseller in your area, please use the contact form for Sales on our homepage www.smartoptics.de.

#### **BRANDS AND TRADEMARKS**

Windows, .NET Framework and Visual C++ are registered trademarks of Microsoft Corporation.

FlyCapture and UpdatorGUI are registered trademarks of Point Grey Research, Inc.

exocad is a registered trademark of exocad GmbH, D-64293 Darmstadt.

Adesso Split and Artist/arTO are registered trademarks of K. Baumann, D-75210 Keltern.

SAM is a registered trademark of SAM Präzisionstechnik GmbH, D-82131 Gauting.

Artex is a registered trademark of Amann Girrbach AG, A-6842 Koblach.

Gamma is a registered trademark of Gamma medizinisch-wissenschaftliche Fortbildungs-GmbH, A-3400-Klosterneuburg.

PROTAR is a registered trademark of KaVo Dental GmbH, D-88400 Biberach/Riß.

Denar is a registered trademark of Whip Mix Corporation, Louisville, KY 40209, USA.

Triple Tray is a registered trademark of Premier Dental, Plymouth Meeting, PA 19462 U.S.A.

Where products belonging to other owners are referenced, we wish to point out that these are brands of the relevant owners that are protected by copyright laws.

Royalty-free icons and illustrations from the Flaticon.com, Fotolia.com, Freepik.com and Pixabay.com image agencies were used in this operating manual.





# Table of Contents

1.	Abou	About this manual		
2.	For yo	For your safety11		
	2.1.	Proper use1	1	
	2.2.	General safety instructions1	2	
	2.3.	Symbols on the device	22	
3.	Scope	e of delivery2	24	
4.	Devic	e components2	29	
	4.1.	Front view	29	
	4.2.	Internal view	31	
	4.3.	Views of the accessories (standard scope of delivery)3	33	
	4.4.	Connections and type plate (rear view)	35	
5.	Setup	and commissioning of the scanner	36	
	5.1.	Choosing a workstation	36	
	5.2.	Unpacking the scanner	36	
	5.3.	Connecting the scanner	11	
	5.4.	Opening and closing the scanner	13	
6.	The object holder system44			
	6.1.	Overview	14	
	6.2.	Matting the surface with 3D scan spray	<b>1</b> 5	
7.	Positi	oning of a single jaw model	16	
	7.1.	Inserting a jaw model with a Multisplit mounting plate4	16	
	7.2.	Inserting a jaw model with an adapter plate	18	
	7.3.	Inserting a jaw model with an object holder	51	
	7.4.	Fastening gingiva, squeeze bites, scanbodies or wax-ups5	54	
8.	Positi	oning of non-articulated occlusion models5	55	
	8.1.	Fastening an occlusion model with a rubber band5	55	
	8.2.	Fastening of an occlusion model with an occlusion clamp. 5	57	

Page 3





9.	Positio	oning of articulated occlusal models	59
	9.1.	Positioning of an occlusal model with an articulator	59
10.	Positio	oning tooth stump models with the multiDie adapter	62
11.	Positio	oning of aTriple Tray® impression	65
12.	Positio	oning objects (universal mode)	69
	12.1.	Fastening objects on the flexible object holder	70
	12.2.	Inserting the object holder	72
	12.3.	Removing the object holder	73
13.	Positio	oning the calibration model	74
14.	Operating principle of the scanner		
	14.1.	Operation via touchscreen	77
	14.2.	Switching between LR and HR mode	79
15.	Device	e maintenance	79
16.	Faults	and repairs	81
17.	Enviro	nmentally-friendly disposal	82
	17.1.	Disposal of the packaging	82
	17.2.	Disposal of the device	82
18.	Techn	ical specifications	83
19.	Declaration of CE conformity8		86
20.	Short	guide to commissioning	87
21.	Glossa	ary	89
22	Indov		Ω1

Document Version: 08-2023 Table of Contents Page





## About this manual

Please read through this operating manual carefully. Keep the operating manual in a safe place. If you pass on the product to somebody else, please make sure that you hand them this operating manual as well.

smart optics accepts no liability for damages due to failure to comply with the operating manual.

#### Introduction

The scanners of the "Vinyl" product family are 3D desktop scanners primarily used to digitize jaw models. Dental laboratories benefit from the precision, speed, durability and great functionality of these devices. Their outstanding technical features include: a fully automatic Z-axis, a 180degree opening that allows open scanning and a touchscreen for easy operation.

Thanks to the large Vinyl system plate, articulators can be scanned; for certain makes condyle-related, otherwise average.

The fully open Vinyl Open Air model provides a solid basic configuration with the most important modules, and it can be expanded when necessary.

The Vinyl standard model is equipped with all modules and also with a lid, allowing closed scanning when necessary.

The high-end model Vinyl High Resolution is the only scanner in the series with a high-resolution camera, a blue light LED and a LED status indicator.

Vinyl Open Air

Vinyl

Vinyl High Resolution

#### **Universal mode**

All scanners of the Vinyl series can be operated in Universal Mode and can therefore be used independently of the dental technology for the highprecision scanning of objects. In the Universal Mode, solid objects that are dimensionally stable, such as tools, jewelry, or toys that are with a maximum size of  $80 \times 60 \times 85$  mm (XYZ) and a maximum weight of 0.938 kg can be scanned.

You can edit the scan data in the open PLY or STL format with most common CAD systems and print the shapes with a 3D printer.

To use the universal mode, dental Scan 3.1.0 or higher is required.





Users of the Universal Mode can limit the reading of this manual to chapters 1 - 6 and 12 - 20.

#### What operating manuals are available?

Generally there are several manuals for smart optics scanners:



#### **Operating manual**

This manual contains the description of the device components, the general notes on safety and the technical data. They will tell you how to operate the scanner and how to handle the accessories.



#### Installation manual

Here you will find out how to install the dental Scan scanner software for the first time and how to receive and install updates at a later date.



#### dental Scan user manual

All scanning methods and the use of optional modules are described in this manual.

#### Graphic symbols and highlighting of important information

In this manual, graphic symbols and text highlights are used to present the information more clearly. In this way you may understand the instructions better.

- A triangle on the left indicates a step involving an action.
- ✓ A tick on the left tells you what the outcome should be if you have followed one or more steps correctly.
- 1. An element in a sorted list
- An element in an unsorted list
- ☐ A list of recommendations or action alternatives

#### Important statements are highlighted in bold type.

Important **KEYWORDS** are highlighted with capital letters in bold type.



Cross-references within a manual are highlighted with a book icon. Crossreferences act as hyperlinks in the electronic version. Click on the number to follow the link. The page navigation of the PDF reader returns you to the starting point ("Previous view" or "Back").

Page 6 Document Version: 08-2023 About this manual









Cross-references to other manuals are highlighted with a symbol for the type of manual: a scanner for the operating manual, a software package for the installation manual and a user on the PC for the dental Scan user manual.



The signpost shows a selection of topics related to your current issue.



This icon points towards functions or modules that are not included in the standard scope of delivery and that can be purchased separately.



Variants in operation and functionality which do not apply to this scanner type are marked with this symbol.



Standard settings specified by the manufacturer are indicated with a steering wheel.



#### TIP

The word "Tip" and the information icon are used to highlight tips on operation and information that require your special attention.



The manufacturer of the product is indicated by the factory symbol.

#### **Figures**

Figures showing equipment and screenshots from the software are used for illustration purposes.

The images shown in this manual may differ from the appearance of your actual device or from the software supplied with your device. The screenshots show realistic examples but do not contain actual data.







#### **Software controls**

The software controls are highlighted as follows:

OK Names of buttons in the software are shown with a colored

background.

**Name** Names of fields in the software are highlighted in bold type.

File | New Names of menus and menu options are separated by a

vertical line.

Name.txt File and path names are shown in Courier font.

#### Units of measurement and numbers

In this manual the numbers are given according to the decimal system with decimal places.

Units of measurement comply with the legal units in Germany as well as the international system of units (SI).

Further information as well as measurement converters can be found on the Internet.

Measurement	Unit	Abbreviation
Mass/weight	Kilogram Gram	kg g
Length	Meter Centimeter Millimeter Micrometer	m cm mm μ
Angle	Degree	0
Time	Hour Minute Second	h m s
Electrical power	Watt	W
Electrical voltage	Volt Alternating current	V AC





Measurement	Unit	Abbreviation
Electrical frequency	Hertz Gigahertz	Hz Ghz
Electrical current	Ampere	А
Protection class	International Protection	IP 1st digit: solids protection; 2nd digit: fluids protection
Temperature	Degree Celsius	0
Memory (PC)	Random-Access Memory	RAM
Data volume	Megabyte Gigabyte	MB GB
CAD data format	Polygon file format Standard Triangulation Language	PLY STL
	American Standard Code for Information Interchange	ASCII e
	Mesh Point Cloud Measurement	MSH PCM
Picture element	Pixel	Px (monitor) dpi (print)
	Megapixel	Мрх





#### Warnings

Warnings provide information on how damage to objects and injury to persons can occur and give instructions on how to avoid risks. Warnings are categorized into four levels depending on the severity of the possible consequences.

#### **NOTE**



This combination of signal word and symbol warns you about possible *material damage* that might occur if the instructions are not followed correctly.

#### **CAUTION**



This combination of signal word and symbol warns you about possible *minor injuries* that might occur if the instructions are not followed correctly.

#### **WARNING**



This combination of signal word and symbol warns you about *severe to lethal injuries* that might occur if the instructions are not followed correctly.

#### **DANGER**



This combination of signal word and symbol warns you of dangerous situations which could lead directly to death or severe injuries.







# For your safety

#### 2.1. Proper use

#### Use in dental technology

The proper use of the "smart optics Vinyl Series" scanners is the optical three-dimensional measurement of human jaw models. The scanners can be used in the field of orthodontics and prosthetics for all kinds of reconstructions, as well as for archiving purposes.

Jaw models in occlusal relationship can be scanned with reference to cephalic presentation, as well as dental registrations (bite registrations) and dental casts (wax-up), and also reference bodies (scan bodies) screwretained in the model.

#### Use as a universal scanner

The proper use of the "smart optics Vinyl Series" scanners is the optical three-dimensional measurement of objects with a maximum size of  $80 \times 60 \times 85$  mm (XYZ). The maximum weight of the objects must not exceed 0.938 kg. The scanners can be used for the production of, for example, tools, jewelry or toys as well as to test the accuracy of replicas.

#### **Material characteristics**

The scanners of the Vinyl series scan materials with a dry, opaque surface of white, saffron, gold, blue, beige, yellow and pink color. Reflective or dark surfaces can be scanned if they are previously treated with 3D scan spray. Moveable parts of the models must be able to be immobilized, so that they do not change their position during scanning.

Document Version: 08-2023 For your safety





#### Improper use

The scanners of the Vinyl series are not designed to scan models made of transparent materials, models of non-dimensionally stable or unstable materials or living organisms.

The scanners of the Vinyl series are not suitable for use in an environment heavily burdened by emissions (e.g. dust or paint). Their use in a private environment is also excluded.

smart optics strongly advises to avoid any other than the proper use. smart optics assumes no liability for damages due to an improper use of the scanner and/or non-observance of the safety instructions indicated in these operating instructions.

## General safety instructions

#### **User qualification**

Users who set up, commission or operate a smart optics scanner require special knowledge for the safe operation of the scanner. This knowledge can be acquired by means of the following measures:

- ☐ Read and follow these operating instructions, in particular the instructions concerning setup, commissioning and cleaning of the scanner.
- ☐ Participate in an induction or training course.
- ☐ Comply with local laws, regulations, occupational safety and accident prevention rules concerning your scanner's workplace.
- ☐ Make sure that all devices and cables pass the periodic safety checks for electrical appliances. Replace damaged devices and cables immediately. Make sure that the technical specifications of power cables are adequate.

#### **Constructive protective measures**

The scanner has been developed and manufactured in compliance with the applicable safety standards and with the greatest possible care to ensure safe operation and to protect the user against injuries.

A fuse to protect the device from overvoltage is integrated into the device.

Document Version: 08-2023 For your safety







#### **Protection against injuries**

Despite the protective measures included in the design, some residual risks that might lead to injuries cannot be excluded. In this section you will learn what measures you can take to protect yourself and other people.

Electric shock and short circuit



#### **WARNING**

#### Risk of injury due to electric shock Fire hazard due to short circuit

A technical defect of the cable or of individual components of the scanner can cause an electric shock or a short circuit. This can result in a fire.

- ▶ Make sure that no electrical equipment can come into contact with water/moisture. Should this however happen at some time, disconnect the power plug immediately. Dry the affected parts with a soft microfiber cloth.
- ▶ Do not work with defective devices or cables under any circumstances.
- Operate electrical equipment only at the recommended operating temperatures.
- ▶ Only use the supplied cable or the corresponding spare parts.
- If electrical devices are not used over a longer period, such as overnight, turn them off and disconnect the power plug from the socket.

Magnetic fields



#### **WARNING**

#### Health hazard due to magnetic fields

The scanner and its accessories include magnetic components. Magnetic fields can be health-damaging.

Persons with implants, in particular heart pacemakers, may only operate the scanner and its accessories with express medical permission.

Document Version: 08-2023 For your safety Page





Stripe light



#### WARNING

#### Health hazard due to stripe light

The scanner works with stripe light. Constant visual contact with stripe light can trigger epileptic seizures, migraine or similar conditions.

Individuals suffering from such conditions should close or cover the scanner during operation.

Blue light



#### **WARNING**

#### Health hazard due to LED light with high blue content

Recent studies show that light with a high blue content can lead to permanent retinal damage or contribute to it. White LED light also has a high blue content. Users should therefore avoid prolonged visual contact with white or blue LED light.

- The scanner emits white light or blue light with each measurement. Close or cover the scanner during operation.
- Wearers of glasses can protect their eyes with glasses with a blue light filter.

Falling



#### **CAUTION**

#### Risk of falling due to packaging materials

The scanner is extensively packaged for protection against transport damages. The packaging can be an obstacle during setup and can cause falls.

- Do not leave packaging materials lying on the ground.
- Before transporting, move all obstacles out of the way.

Document Version: 08-2023 For your safety





Carrying



#### **CAUTION**

#### Risk of injury due to incorrect carrying

Due to its size and weight, we recommend that two people should unpack and set up the scanner. In particular, shorter individuals can be injured if they lift or carry the scanner by themselves without any help.

- Two people should lift the scanner out of its packaging.
- Two people should carry the scanner.
- Carry the scanner with the lid closed.
- To transport the scanner, always grasp it by the bottom corners, never by the lid.

Crushing



#### **CAUTION**

# Risk of injury due to the lid or the mechanical parts of the scanner

You can get crushed by the lid or the mechanical parts of the scanner.

- Always hold the lid by its handle to open or close it. If you act sideways and move the lid, you will crush your fingers between the lid and the housing.
- Lift the lid gently. Otherwise it will fall down again, because there is no resistance in the closing mechanism.
- Always open the lid up to the stop.
- Only put your hands inside the scanner when the axes have stopped. Should the axes not stop at the end of a scanning process, turn off the scanner and disconnect the power plug.

Document Version: 08-2023 For your safety Page 15





Entanglement



#### **CAUTION**

#### Injury due to entanglement of clothing, jewelry or hair

Loose clothing, jewelry or long hair can get entangled in the lid or the mechanical parts of the scanner. Due to motion, objects and hair can get entangled inside the scanner. This can lead to injury.

- Do not wear any loose garments such as shawls and ties, or jewelry like long necklaces, at the scanner's workplace.
- Comb your hair into a firm pigtail, for example.
- However, should a piece of clothing, hair or jewelry get entangled among the moving parts of the device, turn the scanner off immediately. Disconnect the power plug before removing the entangled part.

#### **Protection against material damages**

Despite the protective measures included in the design, some residual risks that might lead to material damages (loss of data or damage to devices) cannot be excluded. In this section you will learn what measures can be taken against material damages.

Optics/Electronics



#### **NOTE**

## Damage to optical and electronic components due to touching, humidity, dirt and cleaning

The optical and electronic components in the interior of the scanner are highly sensitive. Any touch, as well as dirt, moisture and cleaning agents can damage them.

- Do not touch any optical and electronic components,
- Protect the optical and electronic components against soiling by closing the lid or covering the scanner.
- Never use the 3D scan spray inside the scanner.
- If the optical or electronic components have to be cleaned, do not do it yourself. Please contact your reseller.

Document Version: 08-2023 For your safety





#### Carrying



#### **NOTE**

#### Damage to the scanner due to incorrect carrying

The lid and the front upper corners are not suitable to be used as carrying handles. If you carry the scanner by the lid, it can come off and the scanner will fall and get damaged. If you reach into the front upper corners, you can damage sensitive components, especially the camera and the sensor.

- Carry the scanner with the lid closed (with the exception of Vinyl Open Air).
- ▶ To transport the scanner, always grasp it by the bottom corners, never by the lid nor the front upper corners.
- ▶ During transport place a protective film over the lid so that it cannot open.

#### Cleaning agents



#### NOTE

#### Damage to surfaces due to unsuitable cleaning agents

Paper towels, coarse cotton cloth, detergents, polishing paste, etc. leave scratches on the sensitive surfaces.

Only use the recommended materials for cleaning.

#### Moisture



#### **NOTE**

#### Material damages due to moisture

Constant contact with moisture can damage sensitive materials.

- Only operate the scanner at a dry workplace.
- Avoid any contact with water/moisture. Should this however happen at some time, disconnect the power plug immediately. Immediately wipe off any water/moisture using a soft microfiber cloth.

Document Version: 08-2023 For your safety Page





Climate



#### **NOTE**

#### Measuring errors due to unsuitable climatic conditions

The scanner is designed solely for use inside dry, closed rooms. Precise measurement results are only achieved by the scanner under suitable climatic conditions. Too much heat causes measurement errors and overheating of the scanner. Overheating can permanently damage the scanner.

- ▶ Only operate the scanner at temperatures of 18°C 30°C.
- Only operate the scanner in low humidity conditions.
- Avoid direct sun exposure at the workplace.
- Reduce cold, heat and high humidity, e.g. through the use of air conditioners or sun protection.

Reflections



#### NOTE

#### Measuring errors due to reflections on the measured object

Strong extraneous light leads to unwanted reflections on the measured object. This affects the accuracy of the measurements.

- Select a workplace away from a window or not affected by high levels of artificial lighting.
- If no other workplace is available, scan with the lid closed. If you are using Vinyl Open Air, cover the opening during the scanning process, e.g. with light-tight foil.

Vibration



#### **NOTE**

#### Measuring errors due to vibrations

The surface on which the scanner is placed must not vibrate. Vibration will cause inaccurate measurement results.

- Place the scanner on a sturdy, stable base able to support at least twice the weight of the scanner, i.e. 2 x 23 kg.
- Stabilize the base by means of a strutted under frame or by attaching it to a stable wall.





Calibration



#### NOTE

# Inaccurate measurements due to failure to calibrate or use of a damaged calibration model

The measuring accuracy of the scanner is only guaranteed if the scanner is calibrated. This operation requires the use of a calibration model and the associated default values.

The calibration model can be mechanically damaged. This can only be tolerated in the border area.

- Perform the calibration of the scanner after commissioning, after transporting it and whenever the software asks you to do so during operation.
- Only start the procedure when the values registered in the software match the values of the calibration model.
- Check whether the calibration model is damaged at a central location.
- Only use the calibration model if it is in perfect condition.

Model fastening



#### NOTE

# Damage to the scanner and the models due to missing or incorrect fastening

Due to the motion of the axes inside the scanner, unsecured or incorrectly secured models will fall down during the scanning process.

- Do not place any objects on top of or inside the scanner.
- Never place unfastened models inside the scanner.
- Only use the included object holder or approved accessories to fasten models.
- When fastening occlusal models with rubber bands, only use stable, unused rubber bands. Thin or porous rubber bands can tear.
- Always fasten models on the flexible object holder and the multiDie adapter using adhesive pads made of modeling clay.
- Only place stump models prepared with metal pins inside the multiDie adapter.
- ▶ However, should a model still fall down, turn off the scanner immediately. Then remove the model or all broken pieces from the interior of the scanner.

Document Version: 08-2023 For your safety Page







Adhesive materials



#### **NOTE**

#### Damages to models due to unsuitable adhesive materials

Models, object holders and the scanner get soiled or are even damaged through the use of adhesive tape, instant glue or similar products.

Only fasten models using adhesive pads made of modeling clay.

Articulators



#### **NOTE**

# Damages to the scanner and models due to incorrectly positioned articulators

Articulators are not to be fastened inside the scanner. This is possible because articulators can be scanned with reduced axis motion.

- Only place an articulator inside the scanner when the software asks you to do so.
- ▶ Always place an articulator with the front side first inside the scanner.
- Always place an articulator with all legs resting on the system plate.
- Always leave the lid open (except in the case of Vinyl Open Air)
- Remove an articulator immediately when the software asks you to do so.
- ▶ However, should an articulator still fall down, turn off the scanner immediately. Then remove the articulator or all broken pieces from the interior of the scanner.

USB cable



#### **NOTE**

#### Data loss due to a too long USB cable

The transfer of data between the scanner and the PC is affected by the length of the USB cable.

▶ Only use the supplied USB cable or an original spare part with the corresponding length (max. 2 meters).





Magnetic fields



#### **NOTE**

#### Data loss due to magnetic fields

The scanner and its accessories include magnetic components. Metal-containing technical devices and data carriers, e.g. credit cards, can be impaired in their function or even permanently damaged through contact with magnets.

▶ Make sure there is an adequate distance between metalcontaining technical devices/data carriers and the magnet.

#### Response in case of a defect

You must not work with a defective system. Damage to the scanner, its accessories or a faulty function have occurred with high probability if:

- parts have been visibly damaged
- one or more of the moving axes continue to rotate or rotate in an uncontrolled way
- the calibration procedure fails
- measuring errors occur despite performing all work steps
- the software displays an error message.
- Follow the instruction in an error message.
- Chapter 5.3
- ▶ Check the connections. If applicable, replace the USB slots.
- ▶ End the software, switch off the scanner and the PC and restart the system.



- ▶ Check whether the software is installed correctly.
- ▶ Check the calibration data.
- Chapter 15 Follow the instruction for device maintenance.
  - ▶ Update your system. Check whether all important Windows updates have been installed.
  - ▶ Ensure that your computer is free of viruses and malware.
  - Allow the execution of dental Scan in your virus protection software.
- Chapter 16
- If none of these measures solve the problem, repeat with another computer to exclude a computer fault.
- If the problem also occurs with the other computer, refer to your reseller to clarify the cause of the fault.

Document Version: 08-2023 For your safety Page 21





## 2.3. Symbols on the device

#### Warnings



#### Warning against hand injuries

This symbol refers to the lid or the closing mechanism of the scanner. If a body part enters the scanner opening, there is a risk of crushing. Always exercise care when reaching into the scanner to place the object holder.



#### **Electrical voltage warning**

This sign warns of electrical voltage inside the scanner. The rear cover of the scanner may only be removed by skilled personnel if the mains plug has been disconnected from the mains.



#### Do not touch

This symbol refers to the calibrated optical system. This must never be touched or cleaned as it could be damaged in the process.



#### Other

#### **CE** label

The CE label on the type plate documents that the European directives that are applicable to the device were complied with at the time of issue of the CE certificate.

#### **EMC** labelling

EN 55011 Class A



#### **Fuse**

Symbol for an overcurrent protection device, installed in a mains connection, suitable for AC voltage 100 - 240 Volt and a frequency of 50/60 Hertz.

Document Version: 08-2023 For your safety Page 22







#### **USB**

Symbol for a USB connection, also for a camera connection.



## **Protective grounding**

Symbol for connection to an external conductor to protect against electric shock in the event of a fault or connection of protective grounding (mass).



#### On/Off switch

Symbol for a switch to turn the device on and off.





# 3. Scope of delivery

The scope of delivery includes the following components:

Vinyl Open Air (90351)	Vinyl (90334)	Vinyl High Resolution (90400)
<ul><li>1 scanner</li><li>Touchscreen</li><li>White light LED</li><li>1.3 Mpx camera</li></ul>	1 scanner  - Touchscreen  - Front lid  - White light LED  - 1.3 Mpx camera	<ul> <li>1 scanner</li> <li>Touchscreen</li> <li>Front lid</li> <li>Blue light LED</li> <li>LED status bar</li> <li>2.8 Mpx camera</li> </ul>
	Accessory case:	Accessory case:
	<ol> <li>object holder with knurled screw</li> <li>object holder, flexible, with large plate</li> <li>multiDie adapter</li> <li>Triple Tray® impression holder</li> <li>packs of adhesive pads, extra-strong*</li> <li>calibration model</li> <li>data carrier</li> </ol>	<ol> <li>object holder with knurled screw</li> <li>object holder, flexible, with large plate</li> <li>multiDie adapter</li> <li>Triple Tray® impression holder</li> <li>packs of adhesive pads, extra-strong*</li> <li>calibration model</li> <li>data carrier</li> </ol>
Accessory box:	Accessory box:	Accessory box:
<ul> <li>object holder with knurled screw**</li> <li>calibration model</li> <li>data carrier</li> <li>power cable</li> <li>USB cable</li> </ul>	<ul><li>1 power cable</li><li>1 USB cable</li></ul>	<ul><li>1 power cable</li><li>1 USB cable</li></ul>

- \* UHU patafix PROPower or a comparable brand product made of removable, reusable adhesive glue pads
- \*\* To be able to use Vinyl Open Air for all project types, the required optional accessories can be purchased additionally.

Document Version: 08-2023 Scope of delivery Page 24





#### Contents of the data carrier

Vinyl	Vinyl High Resolution	Vinyl Open Air
Calibration data	Calibration data	Calibration data
dental Scan software with the	dental Scan software with the	dental Scan software with the
modules:	modules:	modules*:
<ul><li>secondDie</li></ul>	<ul><li>secondDie</li></ul>	<ul><li>secondDie</li></ul>
– multiDie	– multiDie	<ul> <li>Monochrome texture Scan</li> </ul>
<ul><li>multiCase</li></ul>	<ul><li>multiCase</li></ul>	<ul> <li>Color texture Scan</li> </ul>
<ul> <li>Monochrome texture Scan</li> </ul>	<ul> <li>Monochrome texture Scan</li> </ul>	<ul> <li>Universal mode</li> </ul>
<ul> <li>Color texture Scan</li> </ul>	<ul> <li>Color texture Scan</li> </ul>	
<ul> <li>Triple Tray®impression Scan</li> </ul>	<ul> <li>Triple Tray® impression Scan</li> </ul>	
<ul> <li>Universal mode</li> </ul>	<ul> <li>Universal mode</li> </ul>	
Manuals (PDF):	Manuals (PDF):	Manuals (PDF):
<ul> <li>Scanner Operating Manual</li> </ul>	<ul> <li>Scanner operating manual</li> </ul>	<ul> <li>Scanner Operating Manual</li> </ul>
<ul> <li>dental Scan Installation</li> </ul>	<ul> <li>dental Scan Installation</li> </ul>	<ul> <li>dental Scan Installation</li> </ul>
Manual & User Manual	Manual & User Manual	Manual & User Manual

<sup>\*</sup> Additional modules available for purchase (optional accessories).

#### **Packaging**

The packaging consists of the following components:

- pallet base
- bottom insert with foam protection
- accessory box
- bag (protective cover for the scanner)
- foldable wooden box
- foam insert
- foam pad for the front (Vinyl Open Air only)
- accessory case or box (depending on the model)
- inner lid with handling holes and foam protection
- outer lid





#### **Optional accessories**

In the case of Vinyl and Vinyl High Resolution, all components are already included in the standard scope of delivery. The exception being articulators (not an accessory) as well as adapter plates for articulators (chargeable accessory).

The functionality of Vinyl Open Air can be expanded by purchasing accessories.

You can purchase spare parts and chargeable accessories from your reseller.

Item number	Item name	Description
90404	PC Dell Precision 3630 CTO	High-performance computer with preinstalled software
90073	3D calibration model	Replacement or additional item
90082	Object holder with knurled screw	Replacement or additional item
90350	multiDie adapter	Replacement or additional item
90242	Occlusion clamp	Aid for the occlusion of non-articulated jaw models
90352	Triple Tray® Impression Holder	Replacement or additional item
90358	Object holder, flexible, with large plate	Replacement or additional item
90345	Adapter plate for articulated single jaw models suitable for SAM® AXIOSPLIT®	Additional item

Document Version: 08-2023 Scope of delivery Page 26





Item number	Item name	Description
90346	Adapter plate for articulated single jaw models suitable for GAMMA® Reference	Additional item
90370	Adapter plate for articulated single jaw models suitable for Whip Mix Denar®	Additional item
90371	Adapter plate for articulated single jaw models suitable for KaVo PROTAR®	Additional item
144042	Handmade dust cover	Additional item

For safety reasons, please use only the replacement or additional items mentioned here. For parts that are not listed as replacement items, e.g. cables, please replace them with commercially available products with the same technical specifications.



PCs must meet the technical specifications.



Activate the optional modules in the software that may require special accessories for their use. Further information is given in the dental Scan user manual.

Document Version: 08-2023 Scope of delivery Page 27





#### Supported articulator systems

Articulators are available at specialist dental dealers, but cannot be purchased from smart optics as accessories for the scanner.

All conventional articulators can be used for vestibular scans with articulated occlusion models. Condyle-related measurements are possible with the articulators of the following manufacturers:

- AMANN GIRRBACH Artex®
- Baumann Dental Artist/arTO®
- SAM® AXIOSPLIT®
- GAMMA® Reference
- Whip Mix Denar®
- KaVo PROTAR®


Other articulators can be used for mean value measurements.





# 4. Device components

# 4.1. Front view



#### CAPTION

- 1 HR label (only for Vinyl High Resolution)
- 2 Lid (not included for Vinyl Open Air)
- **3** Device name
- 4 Handle (black or silver (matte) depending on the model/series)
- **5** Touchscreen





#### Detailed view of the touchscreen



#### **CAPTION**

1 Starting status (logo): Operation from the PC possible

**2** Hibernation

(screen saver) Operation from the PC possible

**3** Operating status: Start button and process display

Start button: Operation from the device possible


Process display: The device is working, no user action possible

Document Version: 08-2023 Device components Page 30





## Internal view



#### CAPTION


- 1 Optics (camera and 3D sensor)
- Swivel axis (electromotive lateral movement 0 45°) 2 Z-axis (electromotive up and downwards movement 0 – 30 mm)
- Guide rails for the lid 3 (not present in the Vinyl Open Air variant)
- Rotary axis (electromotive rotation up to 315°) System plate with Adesso Multisplit base plate

Page 31 Document Version: 08-2023 Device components





#### Detailed view of the system plate



#### **CAPTION**

- 1 Turntable with non-slip rubber mat for the placement of articulators
- 2 Adesso® Multisplit base plate for mounting object holders, adapter plates and Multisplit mounting plates

#### Detailed view of camera, sensor and LED



#### **CAPTION**

- 1 LED (white light or blue light (only for Vinyl High Resolution))
- Optical system (camera and sensor) (high resolution for Vinyl High Resolution)

Document Version: 08-2023 Device components Page 3





#### Detailed view of the status-LED (only for Vinyl High Resolution)



**CAPTION** 

Status LED on the swivel axis

■ Purple Hibernation (equivalent to the logo or screen saver on

**1** the touchscreen)

■ Green Operation (corresponds to the start button and the

process display on the touchscreen)

# 4.3. Views of the accessories (standard scope of delivery)



Accessory case with foam insert







Object holder with knurled screw



Flexible object holder with large plate, fixing with adhesive pads





adhesive pads

multiDie adapter, fixing with Triple Tray® Impression Holder



Adhesive glue pads, extrastrong (Example image)



**Calibration model** (Example image)



Power cable



**USB** cable



**Data carrier (Example image)** 

Page 34 Document Version: 08-2023 Device components





# 4.4. Connections and type plate (rear view)



#### **CAPTION**

- 1 ON/OFF switch
- **2** Functional grounding
- 3 USB connection for device control and camera
- 4 Mains connection with fuse

Type plate



3D sensor serial No. SO-20254.00-17-022

The type plate with technical information such as serial number of the device, place and date of manufacture is located on the rear.

The serial number of the 3D sensor is given on a separate label.

Document Version: 08-2023 Device components Page 35





# 5. Setup and commissioning of the scanner

Chapter 20

You will find a short guide to this chapter in the appendix "Short Guide to Commissioning"

## 5.1. Choosing a workstation

Observe the safety instructions!

Chapter 2

- When choosing a workstation, please note that the Vinyl Series scanners are EMC Class A devices suitable for use in a commercial environment.
- Chapter 18
- ▶ Choose a stable work desk that is big enough to allow the scanner to be connected to a PC. The weight of the scanner is 23 kg and its setup dimensions are  $455 \times 430 \times 420$  mm or, with its lid open, approximately  $455 \times 660 \times 470$  mm (W x H x D).
- ▶ Ensure easy access to the rear as this is where the ON/OFF switch is located.





▶ Ensure that there are enough and sufficiently fused mains sockets for all the devices. After consultation with an electrician, you can use an approved extension cable, a multiple power socket and a socket adapter (additional electrical devices not included in the scope of delivery).

## 5.2. Unpacking the scanner

Observe the safety instructions!

Chapter 2

The scanner is delivered in a foldable wooden box. The wooden box is equipped with an impact and tilt indicator (outside) and clad with a foam transport protection (inside).



#### TIP

Keep the packaging for possible transport or return. The original packaging is precisely adapted to the scanner and thus offers the best protection against transport damage. If necessary, you can request packing instructions from smart optics.

Check the outer packaging for visible damage immediately on receipt.





- If an indicator indicates a careless storage or even any shipping irregularities or if damage is evident on the packaging, please inform the delivery company and your reseller. Follow the instructions provided by the delivery company and reseller before continuing to unpack the device.
- Transport the scanner in the packaging as close to the workstation as possible.



To unpack, you needed a flat-tip screwdriver and combination pliers.

- Bend back the metal tabs at the top and bottom of the wooden box. For this purpose, carefully press the screwdriver between the tab and the wood until a gap is created.
- Grasp the tab with the pliers.
- Bend the tab about 180° (those at the top upwards; those at the bottom downwards).
- ✓ When all metal tabs are in a vertical position, they will fit through the slots at the top and bottom.



Unpacking instructions



- Lift the lid and the inside covering.
- Remove the accessory case or box and place it aside. Depending on the model/series, the accessories may be on top or below the scanner.
- Remove the foam board.







- ▶ Pull the wooden box upwards.
- Press the wooden box together from two opposite corners.



### **Carrying by two persons (recommended)**



### TIP

Due to the dimensions and weight smart optics recommends carrying by two persons.

Never hold it by the lid.

In the case of Vinyl Open Air, never hold it by the front upper corners because sensitive technical components are placed there.

The protective foil prevents incorrect gripping or unintentional opening of the lid.

Page 38







- ▶ Position yourself with a second person to the right and left of the scanner.
- ▶ Place one hand each under the front and rear corner of the scanner.
- Lift the scanner simultaneously.
- When carrying, incline the scanner slightly backwards.
- ✓ In this position the lid remains secured.

### Carrying by one person (not recommended)





- Position yourself behind the scanner.
- Grasp with one hand on the side and as far forward as possible under the scanner.
- Lift the scanner so that it is slightly inclined backwards when carrying.
- ✓ In this position the lid remains secured.





### Setting the scanner down

- Set the scanner down carefully and in a straight position at the place of installation.
- ▶ Remove the protective foil.
- If necessary, take the accessory box out of the base insert.



Chapter 3

- Immediately check the scanner and the accessories for completeness and any visible damage.
- ✓ If the scope of delivery is not complete or if one of the supplied parts shows visible damage, inform your reseller immediately. Follow the instructions provided by the reseller before continuing with the setup.





# 5.3. Connecting the scanner

Observe the safety instructions!

Chapter 2

Chapter 4.4

The scanner requires a mains connection and a PC connection. The scanner does not work without a software.

The connections are located on the rear of the scanner.

- Connect the rectangular type B plug of the USB cable to the USB port of the scanner.
- ▶ Connect the device plug of the power cable to the corresponding mains socket on the scanner.



# THE SCANNER

**CABLE CONNECTIONS OF** 

- 1 USB connection
- 2 Mains connection

Connect the mains plug of the power cable to the mains. It is not necessary to use an adapter to adjust the voltage as an AC-DC device is integrated in the scanner.



Connect the flat type A plug of the USB cable to a free USB port of the PC.








- The scanner is operational.
- Switch on the scanner using the power switch.
- ✓ The power switch is in position I.
- ✓ The touchscreen displays the starting status (logo).
- Install dental Scan now. You will find comprehensive information in the dental Scan installation manual.









# 5.4. Opening and closing the scanner


Observe the safety instructions!

Chapter 2

The big lid of Vinyl/Vinyl High Resolution can remain open during operation.

smart optics recommends closing the lid whenever possible (especially during periods of inactivity) to protect sensitive components. If you are using the Vinyl Open Air model, you can cover the opening when necessary, e.g. with light-tight foil.

- ▶ Hold the handle at the lower edge of the lid.
- Lift the lid gently up to the stop.
- Release the lid in this position.



- ✓ The lid will stay in a vertical position without being fastened.
- ▶ To close the lid, grip the handle.
- ▶ Gently move the lid to the end position

Page 43






# The object holder system

#### 6.1. Overview

Observe the safety instructions!

Chapter 2

The extensive object holder system of the Vinyl series allows you to securely fasten and insert jaw models, partial models, tooth stump models, Triple Tray® impressions and occlusion models.



## **CAPTION**


- Object holder: in this case the model with knurled screw for single 1 jaw models
- Base plate with non-slip rubber mat and Multisplit base for object 2 holder, adapter plates and Multisplit mounting plates

Document Version: 08-2023 The object holder system





### Multisplit base plate



#### CAPTION

- 1 Contact points for object holders and adapter plates ("Ovals")
- **2** Contact points for Multisplit mounting plates ("Corners")
- 3 Magnetic adhesive disk for mounting object holders, adapter plates and Multisplit mounting plates

# 6.2. Matting the surface with 3D scan spray

Reflective or dark surfaces must be matted with 3D scan spray before being scanned. smart optics recommends that only products with extrafine grain size and suitability for dental technology be used, for example the scan spray offered as part of the accessories.

- ▶ Follow the instructions of the manufacturer regarding application, safety and material compatibility.
- ▶ If the 3D scan spray has to be sprayed in the work area, close the scanner or cover it previously. This way, you protects the sensitive optical and the electronic components from damaging dusts. Never spray into the inside of the scanner.
- ▶ Fasten the model or object to the object holder before spraying.
- Spray according to the instructions as targeted as possible onto the model or object.

Document Version: 08-2023 The object holder system Page 4





- ▶ When inserting, make sure not to touch the spray layer.
- After scanning, remove the spray residue from the object holder and model or object. Follow the instructions of the manufacturer for this as well.

# 7. Positioning of a single jaw model

You have several options for fastening single jaw models.



- Inserting a jaw model with a Multisplit mounting plate
- Inserting a jaw model with an adapter plate
- Inserting a jaw model with an object holder

# 7.1. Inserting a jaw model with a Multisplit mounting plate

Observe the safety instructions!

Chapter 2

No additional fastening is required for jaw models plastered on a Multisplit mounting plate.

You can place the Multisplit mounting plate directly on the Multisplit base plate of the scanner. This option exists for the articulator systems AMANN GIRRBACH Artex® and Baumann Dental Artist/arTO®.

Plaster the jaw model onto the Multisplit mounting plate according to the manufacturer's instructions. Make sure to insert a magnetic adhesive disk into the Multisplit mounting plate.




▶ Place the Multisplit mounting plate on the Multisplit base plate of the scanner. Ensure that the protrusions ("Corners") on the underside of

<sub>je</sub> 46





the Multisplit mounting plate fit into the recesses of the Multisplit base plate.



- ✓ The Multisplit mounting plate will adhere to the magnet of the Multisplit base plate.
- ✓ The Multisplit mounting plate fits correctly when the plates are exactly congruent. Additional securing of the plastered jaw model is not necessary.



Document Version: 08-2023 Positioning of a single jaw model Page 47





# 7.2. Inserting a jaw model with an adapter plate

### Overview

Using an adapter plate as spacer, you can place the mounting plates of articulated jaw models on the Multisplit base plate of the scanner.



Adapter plates are available as optional accessories for the following articulator systems:









- SAM® AXIOSPLIT®
- GAMMA® Reference
- Whip Mix Denar®
- KaVo PROTAR®

Chapter 7.1

No adapter plates are required for AMANN GIRRBACH Artex® and Baumann Dental Artist/arTO®.


#### **CAPTION**

- 1 Magnet
- 2 Recess ("Corner")
- 3 Protrusion ("Oval")

# Top (in this case: SAM® AXIOSPLIT®)



# Bottom (in this case: SAM® AXIOSPLIT®)



The mounting plates of the corresponding articulator fit the top of the adapter plate. The underside fits to the Multisplit base plate inside the scanner.





### Placing a mounting plate on the adapter plate

Plaster the jaw model onto the mounting plate according to the manufacturer's instructions. Make sure to insert a magnetic adhesive disk into the mounting plate.

Observe the safety instructions! Chapter 2

Place the mounting plate on the adapter plate such, that the protrusions ("Corners") on the underside of the mounting plate fit into the identically shaped recesses of the adapter plate.



- ✓ The mounting plate will adhere to the magnet of the adapter plate.
- ✓ The mounting plate fits correctly when the plates coincide exactly.

SAM® AXIOSPLIT®

► The magnetic adhesive disk in the adapter plates for SAM® AXIOSPLIT® is retained with a fixating screw. Ensure that this is tightened firmly.



✓ Additional securing of the plastered jaw model is not necessary.

### Inserting an adapter plate

Observe the safety instructions!


Chapter 2

- ▶ Hold the adapter plate from the side of the Adesso Split® plastic plate.
- Insert the adapter plate into the scanner.
- ✓ The protrusions ("Ovals") on the underside of the plastic plate engage with the recesses of the Multisplit base plate.

Page 49








- ✓ The adapter plate will adhere to the magnet of the Multisplit base plate.
- ▶ Check whether the adapter plate can be shifted easily. If this is the case, correct the fit until the adapter plate fits securely.

# Removing the adapter plate

▶ Hold the adapter plate on two sides; if necessary, use both hands (not the mounting plate).



- Carefully pull the adapter plate upwards. A certain amount of force is required due to the magnetic attraction.
- ✓ The adapter plate is released from the Multisplit base plate.

Page 50







Service position

### TIP

The axes can be moved unintentionally during removal. For this situation, dental Scan features a function to drive the axes back into the service position.


# 7.3. Inserting a jaw model with an object holder

## Overview

### Object holder with knurled screw

### **CAPTION**

- 1 Rear
- 2 Front
- 3 Floating stop
- 4 Hook
- 5 Threaded pins
- 6 Knurled screw with clockwise thread
- 7 Adesso Split® plastic plate with adhesive disk




### **CAPTION**

- 1 Front
- 2 Rear
- 3 Тор
- 4 Magnetic points
- 5 Adesso Split® plastic plate with adhesive disk

### Flexible object holder

Positioning of a single jaw model







The flexible object holder for Vinyl Open Air is available as an optional accessory from your reseller. It is included in the scope of delivery for Vinyl and Vinyl High Resolution.





### Fastening a jaw model onto an object holder

The object holders for the scanners of the Vinyl series serve to fasten a jaw model mechanically. This method is to be applied to jaw models which are not plastered.

### Object holder with knurled screw

Different sizes of jaw models can be fastened securely to the object holder with the knurled screw.

Observe the safety instructions!

Chapter 2

- Place the jaw model (upper or lower jaw) with the bottom side on the object holder.
- If required, loosen the knurled screw to increase the space.
- ✓ The anterior teeth point in direction of the knurled screw.
- Press the jaw model gently against the threaded pins.
- Tighten the knurled screw.
- ✓ The jaw model is placed correctly if it is flush with the floating stop and the threaded pins.



### Flexible object holder

Jaw models (upper or lower jaw) can be securely fastened on the flexible object holder with the big plate. To fasten them, you require adhesive glue pads.



Two packs of adhesive pads are included in the scope of delivery. You can purchase replacements from office suppliers if required.

- ▶ Cover the topside of the plate with adhesive pads. You should use at least three pads for complete jaw models.
- Place the jaw model (upper or lower jaw) with the bottom side on the adhesive pads.
- ✓ The anterior teeth point in direction of the single magnetic point (1).







- Press the jaw model firmly.
- ✓ The jaw model is placed correctly when it does not extend beyond the magnetic points.
- Tilt the object holder carefully to the right and left.
- ✓ The jaw model is placed correctly when it does not slip.
- ▶ Should the jaw model slip, use additional adhesive pads.


## Inserting the object holders

The described procedure is the same for both object holders.

Observe the safety instructions!

Chapter 2

- ▶ Hold the object holder from the side.
- ▶ Place the object holder with the front side or the knurled screw facing forwards inside the scanner.
- ✓ The protrusions ("Ovals") on the underside of the object holder engage with the recesses of the Multisplit base plate..
- ✓ The object holder will adhere to the magnet of the Multisplit base plate.







▶ Check whether the object holder can be shifted easily. If this is the case, correct the fit until the object holder fits securely.

### Removing the object holder

The described procedure is the same for both object holders. To fasten a jaw model to the object holder, you should always remove the object holder from the scanner.

- ▶ Hold the adapter plate on both sides; if necessary, use both hands.
- ▶ Carefully pull the object holder upwards. A certain amount of force is required due to the magnetic attraction.
- The object holder is released from the Multisplit base plate.



#### TIP



Upon removal, the rotary axis may be moved unintentionally. For this situation, dental Scan features a function to drive the axes back into the service position.

# 7.4. Fastening gingiva, squeeze bites, scanbodies or wax-ups









- ▶ To scan gingiva, squeeze bites, scanbodies or wax-ups, fasten the small prosthetic parts to the jaw model according to the prosthetics manufacturer's instructions.
- ▶ Tighten the scanbodies firmly.

Chapter 7

- ▶ Fasten the prepared jaw model to an object holder or a mounting plate for your articulator.
- ✓ As a rule, no further fixation is required. Loosely placed parts can be glued slightly in two places with a removable adhesive if required.





# 8. Positioning of non-articulated occlusion models

You have various options for inserting an occlusion model into the scanner. The easiest method for non-articulated occlusion is to fasten the occlusion model to the object holder with knurled screw using a rubber band.

The flexible object holder is not suitable for this purpose.

Chapter 8.2

Alternatively you can use an occlusion clamp.

Chapter 9

For articulated occlusions you require an articulator.

# 8.1. Fastening an occlusion model with a rubber band

Observe the safety instructions!

Chapter 2

To fasten the upper and lower jaw in occlusion, you require a conventional rubber band of approx. 0.4 cm width and approx. 8.5 cm diameter. Alternatively you can use a crossed band. The length and strength of the rubber band varies depending on the jaw model. For this reason, always keep several different rubber bands available.

Rubber bands are not included in the optional scanner accessories. For safety reasons, use only new, robust rubber bands and replace these regularly.

- ▶ Place the upper jaw model on the lower jaw model in occlusion.
- ▶ Place the rubber band in cross form so that the top part of the upper jaw model so that an equally long loop hangs from each side.
- Chapter 7.3
- Guide the end of each rubber loop around the hooks on the side of the object holder.



- ✓ The fastening is adequate if the upper jaw model cannot be tipped or shifted with slight pressure.
- If necessary, shorten the rubber bands by coiling over the hooks several times.






Alternatively, you can use two shorter rubber bands or a crossed band. Depending on the length and strength of the rubber bands used, alternative fastening methods are possible, e.g.:

- Place two shorter rubber bands at an angle over the jaw model.
- Guide one end each at the front and rear around one of the hooks on the side of the object holder.

If you are using a crossed band, you can only secure the jaw models to the object holder in occlusion.



Place the crossed band around the jaw models in occlusion such, that one crossed part is on the top and one on the bottom and that the two partial bands are stretched on the side.



- ▶ Secure the connected jaw models on the object holder.
- Guide the laterally stretched rubber bands around the hooks on the object holder.





# 8.2. Fastening of an occlusion model with an occlusion clamp

Observe the safety instructions!

Chapter 2

With the aid of the smart optics occlusion clamps you can fasten the upper and lower jaw model easily and conveniently in non-articulated occlusion. The occlusion is closed in this manner and remains firmly joined, even during mechanical movement.



The occlusion clamp for non-articulated occlusion is an optional accessory for the scanners of the Vinyl series and can be obtained from your reseller.


#### **CAPTION**

- 1 Locking screw
- 2 Mounting rod
- 3 Pressure fitting with ball knob
- 4 Metal pins
- 5 Plug-in connection

### **Overview**

# Occlusion clamp without object holder

# Occlusion clamp with fastened object holder



### **Fastening**



- Fasten first the lower jaw model to the object holder with the knurled screw.
- ▶ Place the upper jaw model on the lower jaw model in occlusion.
- Loosen the ball knob (3) as much as necessary to allow the pressure fitting to be moved forwards and backwards. Fasten the occlusion clamp to the floating stop of the object holder by pushing the metal pins (4) into the holes of the floating stop on the object holder (5).
- Hold the occlusion model firmly.





- ▶ Push the pressure fitting (3) on the ball knob as central as possible over the upper jaw model.
- ▶ Loosen the locking screw (1) until the fastening rod (2) of the pressure fitting can be moved up and down.
- ▶ Hold the ball knob and press the upper jaw model onto the lower jaw model.
- ► Carefully tighten the screw (1) and the ball knob (3) until the occlusion model can no longer be moved.
- ✓ The occlusion model is clamped.



Now place the object holder inside the scanner as usual.





# 9. Positioning of articulated occlusal models

With the help of an articulator, you can establish the articulated occlusion of the upper and the lower jaw model.

You can place any articulator without any further accessories inside the scanner for vestibular scans.

Condyle-related measuring is possible with the following articulators:


- AMANN GIRRBACH Artex® (without adapter plate)
- Baumann Dental Artist/arTO® (without adapter plate)
- SAM® AXIOSPLIT®
- GAMMA® Reference
- Whip Mix Denar®
- KaVo PROTAR®

Articulators are available at specialist dealers, but not from smart optics as an accessory for the scanner.

# 9.1. Positioning of an occlusal model with an articulator

#### CAPTION

- 1 Front
- 2 Rear
- 3 Top support pin
- 4 Front support pin
- 5 Feet







### Inserting the Articulator

# Example

AMANN GIRRBACH Artex®

All articulators are to be handled in the same way, independently of type and manufacturer.



- Make sure that the occlusal model is correctly articulated.
- Remove the top support pin from the articulator.
- Also remove the front support pin, if possible.

Observe the safety instructions!

Chapter 2

It is essential that you wait until the software asks you to insert the articulator. Only then is it ensured that the motion of the axes is minimal during the ensuing scanning procedure.



- ✓ The system plate of the scanner is driven into the horizontal service position. If this is not the case, dental Scan features a function allowing to drive to this position.
- Grasp the articulator by the rear struts with both hands.
- Place an articulator with the front side forward inside the scanner.
- ✓ The front side of the articulator points to the rear side of the scanner. A different orientation is not permitted.

**FALSE:** 

**CORRECT:** 





✓ The articulator stands free and does not touch the scanner.





- Place the articulator on the system plate so that **ALL** feet are resting on the rubber mat. Place the articulator in a central position.
- ✓ The articulator stands centrally on the system plate. A different position is not permitted.
- ✓ The rubber mat is anti-slip. No further fastening is necessary.



#### TIP

Always leave the lid open during an articulator scan.

### Removing the articulator

- ▶ Wait until the software reports the end of the scanning procedure.
- Remove the articulator as soon as the software asks you to do so.
- Grasp the articulator as before by the rear struts.
- Lift the articulator horizontally out of the scanner.
- Now you can proceed with the workflow in the software.



### TIP



Upon removal, the rotary axis may be moved unintentionally. For this situation, dental Scan features a function to drive the axes back into the service position.





# 10. Positioning tooth stump models with the multiDie adapter

You can position up to twelve tooth stump models individually with the multiDie adapter.

You can use the multiDie adapter for the two modules multiDie and multiCase.




The multiDie adapter for Vinyl Open Air is available as an optional accessory from your reseller. It is included in the scope of delivery for Vinyl and Vinyl High Resolution.

### Overview

#### **CAPTION**

- 1 Adesso Split® plastic plate with adhesive disk
- 2 Slots
- 3 Flat side (front)



### Fastening tooth stump models

To ensure that the tooth stump models sit firmly on the multiDie adapter, you require the adhesive pads as fixing material.

Two packs of adhesive pads are included in the scope of delivery. As replacement you can use commercially available adhesive glue pads. These should be extra strong and removable, not coloring nor hardening.

- Fill all slots with adhesive pads.
- ✓ The glue material should fill the slots up to their edge, but nor protrude from them. The glue material can remain in the slots permanently.
- Insert the metal pins of the tooth stump models into the slots.







Release (multiDie)
Project definition multiCase

### TIP

If you only insert of 1, 2 or 3 tooth stump models, the slots 2, 5 and 8 must be used. For all further tooth stump models you can choose the slots freely.

### **Correct fastening**

The prepared tooth stump models are fastened correctly if:

- they are straight (not tilted)
- they have no contact with their neighbor
- they lie on the multiDie adapter in other words, insert the metal pins as deep as possible into the glue material
- they remain in their position when the multiDie adapter is rotated or tilted.

It is of advantage, but not a precondition, to insert the tooth stump models such that the buccal side faces outwards.

### Inserting the Triple Tray® impression holder



Observe the safety instructions!

Chapter 2

- ▶ Hold the multiDie adapter from the side.
- ▶ Place the multiDie adapter inside the scanner so that the flat side faces the front of the scanner.
- ✓ The protrusions ("Ovals") on the underside engage with the recesses of the Multisplit base plate..





- ✓ The Triple Tray® impression holder will adhere to the magnet of the Multisplit base plate.
- ▶ Check whether the multiDie adapter can be shifted easily. If this is the case, correct the fit until the multiDie adapter fits securely.

### Removing the Triple Tray® impression holder

- ▶ Hold the multiDie adapter from the side, if necessary, use both hands.
- ▶ Carefully pull the Triple Tray® impression holder upwards. A certain amount of force is required due to the magnetic attraction.
- ✓ The Triple Tray® impression holder is released from the Multisplit base plate.



### TIP



Upon removal, the rotary axis may be moved unintentionally. For this situation, dental Scan features a function to drive the axes back into the service position.





# 11. Positioning of aTriple Tray® impression

The Triple Tray® impression holder allows you to position Triple Tray® impression trays inside the scanner.

You can use the Triple Tray® impression holder with the Triple Tray® impression scan module.



The Triple Tray® impression holder for Vinyl Open Air is available as an optional accessory from your reseller. It is included in the scope of delivery for Vinyl and Vinyl High Resolution.

### Overview

#### **CAPTION**

- 1 Upper part
- 2 Lower part
- 3 Rear
- 4 Front
- 5 Magnetic adhesion point
- 6 Clamp with spring
- 7 Adesso Split® plastic plate with adhesive disk



### Fastening a Triple Tray ® impression scan

Observe the safety instructions!

Chapter 2

A Triple Tray® impression tray can be clamped quite easily in the Triple Tray® impression holder.

▶ Press the clamp on the broad side and keep it pressed.









- The opposite side is pressed apart.
- ▶ Turn the side of the Triple Tray®impression demanded by the software upwards.
- Insert the flat shaft of the Triple Tray® impression holder into the narrowest position between the open side of the clamp.
- Release the clamp.



### **Correct fastening**

The Triple Tray® impression is fastened correctly if:

- it is clamped
- it is positioned parallel to the lower part of the Triple Tray® impression holder
- the impression to be scanned faces upwards.

### Turning a Triple Tray® impression


To scan the other side of the Triple Tray® impression, turn the upper part without releasing the impression from the clamp.

- Separate the upper part from the lower part applying a little force.
- ✓ The magnetic contact is released.
- ▶ Turn the upper part so that the other side of the impression faces upwards.





Press the lower magnetic point against the lower part.



The upper part will adhere magnetically to the lower part. The adhesive effect is reinforced by the grooves and braces.

### Inserting the Triple Tray® impression holder



- Hold the Triple Tray® impression holder from the rear.
- Place the Triple Tray® impression holder with the impression facing forwards inside the scanner.
- ✓ The Triple Tray® impression holder will adhere to the magnet of the Multisplit base plate.
- ✓ The protrusions ("Ovals") on the underside of the plastic plate engage with the recesses of the Multisplit base plate.
- The Triple Tray® impression faces to the rear of the scanner.

Correct False











▶ Check whether the Triple Tray<sup>®</sup> impression holder can be shifted easily. If this is the case, correct the fit until the Triple Tray<sup>®</sup> impression holder is adhered securely.

### Removing the Triple Tray® impression holder

- ▶ Hold the Triple Tray® impression holder by the lower part; if necessary, use both hands.
- ▶ Carefully pull the Triple Tray® impression holder upwards. A certain amount of force is required due to the magnetic attraction.
- ✓ The Triple Tray® impression holder is released from the Multisplit base plate.



### TIP



Upon removal, the rotary axis may be moved unintentionally. For this situation, dental Scan features a function to drive the axes back into the service position.





# 12. Positioning objects (universal mode)


Scannable objects are very individual physical bodies. Thus, only general advice on correct positioning can be provided here. In principle, smart optics recommends using the flexible object holder for the Universal Mode.

#### **CAPTION**

- 1 Front
- 2 Rear
- 3 Тор
- 4 Magnetic points
- 5 Adesso Split® plastic plate with adhesive disk


### Flexible object holder





With this object holder, solid, opaque objects such as tools, jewelry or toys with a maximum size of  $80 \times 60 \times 85$  mm (XYZ) and a maximum weight of 0.938 kg can be scanned.

 $80 \times 60 \times 85 \text{ mm (XYZ)}$ 0.938 kg





The flexible object holder for Vinyl Open Air is available as an optional accessory from your reseller. It is included in the scope of delivery for Vinyl and Vinyl High Resolution.

Please note that the standard scope of delivery of Vinyl and Vinyl High Resolution includes one multiDie adapter and one Triple Tray® impression holder intended exclusively for dental applications.





# **Properties of scannable objects**

Chapter 6.2

Glossy surfaces (such as metal) or very dark surfaces cannot be scanned. Reflections and lack of contrast to the scanner's interior make an accurate measurement impossible. Therefore, such objects should be matted with 3D scan spray.

### **Dimensional stability**

Only objects that do not change their shape during scanning can be scanned. Make sure that moving parts are adequately immobilized.

### Top and bottom

To fully digitize a three-dimensional object, it usually needs to be scanned from two sides. The second scan captures the bearing surface that was hidden during the first scan. In the simplest case, turn the object on the object holder. Depending on the shape of the object, it may be useful to choose a different fastening for the second side.

# 12.1. Fastening objects on the flexible object holder

Observe the safety instructions!



Objects can be fastened securely on the flexible object holder with the big plate.



Chapter 7.3, 10, 11

#### TIP

Try out the best way to fasten an object.

Test the stability of the fastening outside the scanner using sideways and rotational movements to determine which object holder best fits an object. Create sample scans and compare the results.

So you can, for example, place jewelry instead of single-tooth models into the multiDie adapter or clamp them into the Triple Tray® impression holder. The object holder with the knurled screw can achieve good results with solid bodies having a wide, smooth underside. Please note that multiDie adapters and Triple Tray® impression holders can only support smaller, lighter weight items (<80 x 60 x 85 mm (XYZ), 0.938 kg).

Document Version: 08-2023 Positioning objects (universal mode)





### **Magnetic objects**

Four magnetic points are placed on the top of the object holder. Each of the magnets measures:  $10.0 \varnothing \times 3.0$  H mm. The adhesive force depends on the opposite pole and the environmental conditions and must therefore be checked on a case-by-case basis.

#### Magnetic points



- Place a magnetic object on the plate so that it is held by the magnetic points.
- If the object does not cover the area between the magnetic points, you can enlarge the area with magnetic plates or magnetic tape (available as office supplies).
- ✓ The adhesive effect depends on the individual magnetization of the object and the inserted magnetic plates.

# Non-magnetic objects



Adhesive glue pads are necessary to fasten non-magnetic objects. Two packs of adhesive pads are included in the scope of delivery. You can purchase replacements from office suppliers if required.

- ▶ Cover the topside of the plate with adhesive pads. Different numbers of adhesive pads may be needed, depending on the size, weight and shape of the object. Knead several adhesive pads together and shape e.g. a plate or a pillow.
- ▶ Place the object with the underside on the adhesive pads.
- Align the front / main side of the object with the single magnetic point(1).

Document Version: 08-2023 Positioning objects (universal mode) Page 7







- Press the object firmly.
- ✓ The object is placed properly if it does not protrude beyond the magnetic points.
- Tilt the object holder carefully to the right and left.
- ✓ The object is placed correctly when it does not slip.
- If the object does slip, use more adhesive pads or alternatively test another object holder.

# 12.2. Inserting the object holder

The described procedure is the same for all object holders.


Observe the safety instructions!

Chapter 2

- ▶ Hold the object holder from the side.
- ▶ Place an object holder with the front side forward inside the scanner.
- ✓ The protrusions ("Ovals") on the underside of the object holder engage with the recesses of the Multisplit base plate..
- ✓ The object holder will adhere to the magnet of the Multisplit base plate.







▶ Check whether the object holder can be shifted easily. If this is the case, correct the fit until the object holder fits securely.

## 12.3. Removing the object holder

The described procedure is the same for all object holders. To fasten an object to the object holder, you should always remove the object holder from the scanner.

- ▶ Hold the adapter plate on both sides; if necessary, use both hands.
- ▶ Carefully pull the object holder upwards. A certain amount of force is required due to the magnetic attraction.
- ✓ The object holder is released from the Multisplit base plate.



#### TIP



Upon removal, the rotary axis may be moved unintentionally. For this situation, dental Scan features a function to drive the axes back into the service position.

Document Version: 08-2023 Positioning objects (universal mode) Page 73





## 13. Positioning the calibration model

A calibration model is required to calibrate the scanner. A calibration model is included in the scope of delivery. Each calibration model is measured industrially and assigned these measurements.



#### TIP


If your calibration model is damaged and you are unsure whether this model can be used, your reseller will be pleased to advise you and provide you with a replacement if necessary.

It is recommended to have the calibration model always at hand during operation.

#### **CAPTION**

- 1 Top
- 2 Bottom
- 3 Depressions
- 4 Measured values
- 5 Model number

### Overview








### Inserting and removing the calibration model

Chapter 7.3

The calibration model is to be handled in the same manner as a single jaw model. To fasten it, you require the object holder with the knurled screw.



Observe the safety instructions!

Chapter 2

- If required, loosen the knurled screw to provide more room on the object holder.
- ▶ Place the underside of the calibration model on the object holder.
- ▶ Press the depressions of the calibration model slightly against the threaded pins of the object holder.
- Tighten the knurled screw.
- ✓ The calibration model is placed correctly when it is flush both at the front and the rear.
- ✓ The impression with the measured values faces the floating stop.
- Insert the object holder. When doing it, observe the software's messages.



TIP

Calibration

Information on calibration is provided in the dental Scan user manual.

Page 54

Remove the object holder as usual.

Document Version: 08-2023 Positioning the calibration model Page





## 14. Operating principle of the scanner




#### TIP



Information on all scanning procedures, the setting options and the use of optional modules can be found in the dental Scan user manual.

The most important components of the scanner are the 3D sensor and a positioning mechanism.



The positioning mechanism in the interior of the scanner consist of a rotary and swivel axis driven by an electric motor(4) as well as an automatic z-axis (2).

The freely rotating Multisplit base plate (3) positions the object to be measured with regard to the 3D sensor (1), which is located above the swivel axis.

Document Version: 08-2023 Operating principle of the scanner Page 7







The swivel axis moves the rotary axis with the object holder to the side so that the 3D sensor can capture the object to be measured from the side.

During a measurement, the 3D sensor projects a striped pattern onto the object being scanned.

The light stripes are generated with a white or blue LED light. A blue LED light is used in combination with a 2.8 Mpx camera, as it requires more brightness and sharper contrasts.

The camera records the striped pattern. With the aid of several camera images taken from different angles, dental Scan calculates a 3-dimensional image of the object.

## 14.1. Operation via touchscreen

The scanner is equipped with a touchscreen which allows you to start scanning processes. Operation via touchscreen is an alternative to a mouse click in the software.



You can also control the scanner exclusively via the software. You will find further information in the dental Scan user manual.

The company logo is displayed after switching on the scanner with the software already running.



After starting the software with the scanner already switched on, an animated red circle will be displayed:








▶ When the software displays a prompt, a start button appears on the touchscreen.



- ▶ Follow the instructions given by the software.
- Press the button to start.



- ✓ The scanning process starts.
- ✓ While scanning, the animated red circle will be displayed, and when the procedure has finished, the company logo.
- ✓ After prolonged inactivity, the screen saver is activated automatically.







## 14.2. Switching between LR and HR mode

The software allows you to switch between the low-resolution mode and the high-resolution mode. This determines the level of detail of the subsequent measurements. Please note that very detailed measurements take longer. If this is not required for the current project, you can save time with the low-resolution mode.



The high-resolution mode and thus the switching option is only available for the Vinyl High Resolution.



Whether high or low resolution should be default is determined in the settings of dental Scan. However, you can switch between the modes for each project and every single scan. The exact procedure is described in the dental Scan user manual.

### 15. Device maintenance

The scanner is a delicate and sensitive optical device. To ensure trouble-free operation of the scanner, it is necessary to regularly follow the correct care measures.

Further maintenance measures by the user are not required.

#### **Scanner calibration**

The scanner is calibrated with the aid of the software. This procedure is necessary to ensure accurate measurements.



Perform the calibration of the scanner after the first setup, every time after transporting it and whenever the software asks you to do so during operation. Follow the dental Scan user manual for this purpose.

Document Version: 08-2023 Device maintenance Page 79





### Cleaning the scanner

Observe the safety instructions!

Chapter 2

The scanner should be cleaned regularly during operation. However, never clean the optical or other electronic components.

If the optical or electronic components are soiled, please contact your reseller.

Before cleaning

- Switch off the scanner for safety reasons.
- Remove the cable connections.
- Remove the object holder.

Materials

- Use microfiber cloths which are expressly suitable for high-gloss surfaces. The material must be soft, smooth, lint-free and anti-static.
- Always dampen the cleaning cloth with a glass cleaner.
- Never clean sensitive surfaces with abrasives, polishing pastes or rough cloths.

Cleaning

Document Version: 08-2023

- ▶ Clean the surfaces as well as the object holder, the system plate and the axes with the moistened cleaning cloth.
- ▶ Clean the touchscreen without applying pressure.
- Coarse dirt and debris can be removed from the scanner floor with a vacuum cleaner. Attach the crevice nozzle and set the vacuum cleaner to its lowest vacuuming level.

Device maintenance





## Faults and repairs

Chapter 2, Page 21

If a fault occurs, first follow the safety instructions concerning response to a defect. If the fault persists, contact your reseller to clarify the cause of the fault.

Only have repairs carried out by your reseller.

If you are looking for a reseller in your area, please use the contact form for Sales on our homepage www.smartoptics.de. Improper repairs can cause the scanner to stop working perfectly.

Warranty

smart optics accepts no liability for damages caused by improper repairs. Please note that in this case your warranty claim will also become void.

#### **Serial numbers**

If you have questions or complaints relating to your device, please have the serial number of your Scanner and the serial number of the 3D sensor at hand.

You will find these numbers on the rear of the device:

#### Type plate

Serial number of the scanner



#### Label

Serial number of the 3D sensor

3D sensor serial No. SO-20254.00-17-022





## 17. Environmentally-friendly disposal

The information in this chapter refers to EU directives and German law.

In non-European countries you must follow the corresponding national legislation for the disposal of packaging and electronic scrap.

You can prevent negative consequences for people and avoid harming the environment by the proper disposal of the device.

## 17.1. Disposal of the packaging

In accordance with the EU Directive 94/62/EC | 2015/720/EU and the German Packaging Act (VerpackG), you can return the packaging to your reseller for disposal in Germany. However, smart optics recommends that you keep the packaging in case you need it to transport the scanner or to send it back in the event of warranty claims.

## 17.2. Disposal of the device



The devices marked with this symbol are subject to European Directive 2002/96/EC for WEEE (Waste Electrical and Electronic Equipment).

WEEE registration number of smart optics: DE47893210

Electrical equipment must not be disposed of with domestic waste.

Please note that the scanner is a device that only serves for commercial or industrial use. Disposal via public waste management authorities is therefore not possible.

The device must be returned to the manufacturer for disposal. If you are resident within the area in which the EU directive applies you can also return the device to your reseller.

Document Version: 08-2023 Environmentally-friendly disposal Page 8





# 18. Technical specifications

	Vinyl Open Air	Vinyl	Vinyl High Resolution
FIRST PRODUCTION	2017	2017	2018
LAST HARDWARE UPDATE	2017	2017	2018
HOUSING			
Dimensions W × H × D mm closed open	455 × 430 × 420 455 × 430 × 420	455 × 430 × 420 455 × 660 × 470	455 × 430 × 420 455 × 660 × 470
Weight	23 kg	23 kg	23 kg
Maximum load-bearing capacity (Universal Mode)	0.938 kg	0.938 kg	0.938 kg
Axes	1 Rotary axis 1 Swivel axis 1 Z-axis	1 Rotary axis 1 Swivel axis 1 Z-axis	1 Rotary axis 1 Swivel axis 1 Z-axis
Material	Metal, plastic (ABS HB PMMA)	Metal, plastic (ABS HB PMMA)	Metal, plastic (ABS HB PMMA)
Housing color	black-white (highly glossy)	black-white (highly glossy	)Black-white (highly glossy)
Handle and display frame color	silver matte	silver matte	black matte
Lid	No	Yes	Yes
Touchscreen	Yes	Yes	Yes
ON/OFF switch	Toggle switch, rear	Toggle switch, rear	Toggle switch, rear
Status LED	No	No	Yes

Document Version: 08-2023 Technical specifications Page





	Vinyl Open Air	Vinyl	Vinyl High Resolution
E-TECHNOLOGY			
Supply voltage	100 - 240 V AC 50/60 Hz	100 - 240 V AC 50/60 Hz	100 - 240 V AC 50/60 Hz
Fuse	2 × T 1.6 A L 250 V	2 × T 1.6 A L 250 V	2 × T 1.6 A L 250 V
Power consumption	max. 60 W	max. 60 W	max. 60 W
Connections	1 × USB 1 × power	1 × USB 1 × power	1 × USB 1 × power
Cables & plugs USB	1.8 m   A/B	1.8 m   A/B	1.8 m   A/B
Mains	2.5 m   AC 110/230 V   E+F CEE 7/7	2.5 m   AC 110/230 V   E+F CEE 7/7	2.5 m   AC 110/230 V   E+F CEE 7/7
TEMPERATURE			
Operation:	18°C - 30°C	18°C - 30°C	18°C - 30°C
Storage	-5°C - 50°C	-5°C - 50°C	-5°C - 50°C
MEASUREMENT			
High-resolution camera	No	No	Yes
HR Mode resolution LR Mode resolution Camera	— 1,3 Mpx Point Grey	— 1,3 Mpx Point Grey	2.8 Mpx 1.4 Mpx Point Grey
3D measuring technology	Stripe light triangulation with white light LED	Stripe light triangulation with white light LED	Stripe light triangulation with blue light LED
Measuring field XYZ mm	80 × 60 × 85	80 × 60 × 85	80 × 60 × 85
Measuring accuracy According to ISO 12836	6 µm	6 μm	4 μm
MeasurementScanning of complete Processing jaw Overall time	LR Mode 16 seconds 13 seconds 29 seconds	LR mode 16 seconds 13 seconds 29 seconds	HR Mode* 18 seconds 17 seconds 35 seconds * LR Mode as Vinyl

Document Version: 08-2023 Technical specifications Pag







		Vinyl Open Air	Vinyl	Vinyl High Resolution
Measuremen of single stump	tScanning Processing Overall time	LR mode 33 seconds 08 seconds 41 seconds	LR mode 33 seconds 08 seconds 41 seconds	HR mode* 35 seconds 14 seconds 49 seconds * LR Mode as Vinyl
Measuremen of 3-unit bridge	<sup>t</sup> Scanning Processing Overall time	LR mode 45 seconds 22 seconds 67 seconds	LR mode 45 seconds 22 seconds 67 seconds	HR mode* 50 seconds 25 seconds 75 seconds * LR mode as Vinyl
Measuremen universal Ø	Scanning Processing Overall time	45 seconds 22 seconds 67 seconds	45 seconds 22 seconds 67 seconds	50 seconds 25 seconds 75 seconds

## SYSTEM REQUIREMENTS (PC & CAD-SOFTWARE AVAILABLE EXTRA)

PC (Minimum)	CPU: i3 with 4 × 3.6 GHz RAM: 16 GB Port: USB 3.0 HDD: 40 – 100 GB Graphics card RAM: 1 GB	CPU: i3 with 4 × 3.6 GHz RAM: 16 GB Port: USB 3.0 HDD: 40 – 100 GB Graphics card: RAM 1 GB	CPU: i5 with 4 × 3.8 GHz RAM: 32 GB Port: USB 3.0 HDD: 100 – 250 GB Graphics card RAM: 6 MB
PC (recommended)	CPU: i7 with 6 × 4.7 GHz RAM: 32 GB Port: USB 3.0 HDD: 100 -250 GB Graphics card RAM: 6 MB	CPU: i7 with 6 × 4.7 GHz RAM: 32 GB Port: USB 3.0 HDD: 100 -250 GB Graphics card RAM: 6 MB	CPU: i7 with 6 × 4.7 GHz RAM: 32 GB Port: USB 3.0 HDD: 100 -250 GB Graphics card RAM: 6 MB
Operating system (Minimum)	Windows 7 (64-Bit)	Windows 7 (64-Bit)	Windows 10 (64-Bit)
Operating system (recommended)	Windows 10 (64-Bit)	Windows 10 (64-Bit)	Windows 10 (64-Bit)
dental Scan	as of version 3.0.0	as of version 3.0.0	as of version 3.1.0
CAD-Software	exocad® DentalCAD (all versions)	exocad® DentalCAD (all versions)	exocad® DentalCAD (all versions)

Document Version: 08-2023 Technical specifications Page





## 19. Declaration of CE conformity



smart optics Sensortechnik GmbH Lise-Meitner-Allee 10 44801 Bochum, Germany

## **Declaration of CE- Conformity**

According to EU-regulation 2006/42/EG Appendix II A

We declare that the device identified below complies with the requirements of the EU guideline which regard to safety and physical health requirements both in concept and construction put in circulation.

This declaration becomes invalid in case of an unauthorized change of the device.

Device description: Optical 3D scanner

Device type: Vinyl

EU guidelines applicable: machine guideline (2006/42/EG)

low voltage guideline (2014/35/EU) EMC guideline (2014/30/EG)

Harmonized standards applied:

DIN EN ISO 12100:2010 Safety of machinery

DIN EN 61326-1:2013 DIN EN 61010-1:2010

The CE label was used first for this product in 2017.

Document prepared by: Jörg Friemel

Bochum, 22.05.2017

smart optics
Sensortechnik GmbH
Lise-Meitner-Allee 10
D-44801 Bochum / Germany
Fon: +49 234 29 82 8-0 Fax: -20

Document Version: 08-2023 Declaration of CE conformity Page 86





## 20. Short guide to commissioning

Chapter 5

Detailed work steps and explanations: Setup and Commissioning of the Scanner

#### Workstation





Table's load-bearing capacity: 46 kg Scanner's space requirements  $B \times H \times D$  mm:  $455 \times 430 \times 420 - 455 \times 660 \times 470$ 



PC with recommended configuration

目

CPU: i7 with  $6 \times 4.7$  GHz

RAM: 32 GB Port: USB 3.0

HDD: 100 -250 GB



Graphics card 6 MB RAM

Windows

10 64 Bit

#### Unpacking





▶ Tools required





Observe the display and the instructions



- Follow the description
- Keep the packaging







### Carrying

Step 4











Hold the device under its edges



▶ Tilt the device slightly backwards

### Connection





- Connect the power cable
- Connect the USB cable



Switch on (rear)

## Installation







Copy the data carrier locally





Install the software.



- Start the software.
- Load the calibration data.





- Calibrate: Fasten the model
- Insert #1 and #2



Ready to scan!





# 21. Glossary

3D Scan	Three-dimensional image of the model in the software.
3D sensor	Electronic component for three-dimensional measurement. The 3D sensor is not movable.
Articulator	Dental tool for fabricating a condyle-related occlusion. Articulators are available from various manufacturers. The most common articulators can be scanned condyle-related with the scanners of the Vinyl series.
Blue light/white light	LED light used for measurements.
Calibration	A metrological term.
	A) alignment of the scanner to the values of an industrially calibrated model.
	B) alignment of the motion axes with regard to the individual calibration data.
Calibration model	Industrially measured model which is used for calibrating the scanner.
Height alignment	Positioning of the model in the measuring field of the 3D sensor with the aid of the z-axis.
HR mode	High-resolution scanning.
Lens	Optical component of the camera.
Measurement	Calculation of the surface which can be measured through the projection of stripe light.  Triangulation is the measuring principle.
Measuring field	Size of the maximum area which can be captured by the 3D sensor.
Model	A model which is scanned, e.g. the impression of a plaster jaw.

Document Version: 08-2023 Glossary Page 89





Multisplit base plate	Permanently mounted magnetic plate onto which the Multisplit mounting plates, object holders and adapter plates can be fastened.
Object holder	A holder on which the object to be measured (e.g. jaw model) is mounted and which itself is then mounted to the base in the scanner.
Occlusion clamp	Special holder for smart optics scanners with which a non-articulated jaw model is placed inside the scanner. The occlusion clamp is used for easy fastening of occlusion models.
Rotary axis	One of the motion axes of the scanner. The base can be rotated.  Free rotation of the base allows full circumferential positioning in front of the camera.
Stripe light	A striped pattern which is projected on the model to measure the surface three-dimensionally.
Swivel axis	One of the motion axes of the scanner. The swivel axis moves sidewards so that the model is positioned in front of the camera at different angles.  The swivel axis carries the rotary axis.
Triangulation	Measuring method for determining the position of a point in space with the help of triangles.
Universal Mode	Use of the scanner for non-dental purposes.
Z-axis	One of the motion axes of the scanner. The z-axis moves up and down so that the model is positioned in front of the camera at different heights.

Document Version: 08-2023 Glossary Page 90





## 22. Index

μ		
μ (Unit of measurement)	 	8
3		
3D scan spray Use		
3D sensor		
Α		
Accessories	 	36
Scope of delivery		
Adapter		
ArticulatormultiDie/multiCase		
Scope of delivery		
Triple Tray®		
Adhesive pads		
Application (as intended)		
Articulated jaw model		
Articulator		
Insert/remove		
Manufacturer	,	
Single jaw model		
ASCII (data format)		
Axes	 	31
В		
Blue light	 	76
Box (Transport)	 	36
C		
Cable		/ <sub>1</sub> 1
Cabling		
CAD software		
Calibration	 	79
Calibration data		
Calibration model		
Case (accessories)		
Caution (convention)	 	10
Conformity	 	86
Labelling		

Clamp	
Occlusion	57
Triple Tray®	65
Clasp forceps	57
Cleaning	80
Cleaning agents	80
Code of response	21
Colors	
Status LED	33
Colors (material)	11
Comma	
Notation of numbers	8
Connections	
Mains35,	41
USB35,	41
Conventions	6
Cover	43
D	
Data carrier	24
Defect	
Calibration model	74
Device 21,	81
Denominations	
Accessories	33
Adapter plate	48
Articulator	59
Calibration model	74
Exterior view	29
Internal view 31, 32,	44
multiDie Adapter	62
Object holder	51
Occlusion clamp	57
Rear view	35
System plate	45
Touchscreen	30
Triple Tray® impression holder	65
dental Scan77,	79
Scope of delivery	24
Depth of the device	83
Dimensions (technical data)	83
DIN standards	86
Display 30,	77
Disposal	
Device	82
Environment	82
Packaging	82





E		Fastened with an articulator	59
_		With gingiva	54
Electromagnetic compatibility	36	With scanbody	54
EMC		With squeeze bite	54
Directive	86	With wax-up	54
Labelling	22		
Workstation	36	K	
Environment	12, 82		
Error	21, 81	kg (technical data)	83
EU Directives	86	kg (unit of measurement)	8
		Knurled screw	52, 57, 74
F			
		L	
Fastening the calibration model inside the sc		/	
Fault		LED (white/blue)	
Fuse		Lifting (scanner)	
Fuse (technical data)	84	Lifting points	
		Low Voltage Directive	
G		LR and HR mode	79
Gingiva	54	N/I	
Gingival sculpture		M	
Grounding or protective grounding		Machine Directive	86
drounding or protective grounding	23	Mains connection	
ш		Maintenance	
н		Manuals (types)	
Handle (Lid)	43	Material (colors)	
Height of the device		Material damages (Protective measures)	
High Resolution (Vinyl)		Matting (surface)	
Holder	= .	Measured values (Calibration model)	
Articulator	59	Measuring accuracy	
Mounting plate		Measuring field (technical data)	
multiDie/multiCase		mm (Unit of measurement)	
Object holder		Model holder	
Occlusion		Mounting plate	
Triple Tray®		Connect to adapter plate	49
Hz (technical data)		Direct insertion	
		MSH (data format)	
I .		multiDie/multiCase	
1		Adhesive pads	62
Impression		Fasten tooth stumps	
Clamping	65	Module	
Insert/remove		Scope of delivery	24
Injuries (Protective measures)	13	Slots	
Inserting tooth stumps individually			
Installation		N	
I		Note (convention)	10
-		Number (Calibration model)	74
Jaw model		Numbers	
Fasten with occlusion clamp	57	Decimal system	8
Fasten with rubber band	55		
Fastened inside the scanner 46, 48,	52, 55, 60		





### 0 Object Fastening inside the scanner ...... 70 Object holder...... 51, 69 Articulator......59 Fasten gingiva .......54 Fasten scanbody......54 Fasten squeeze bite ...... 54 Fasten wax-up......54 Occlusion clamp ...... 57 Rubber band ...... 55 Objects (Universal Mode) ...... 70 Occlusion......59 Fasten the model with an articulator...... 59 With occlusion clamp ...... 57 With rubber band......55 Occlusion clamp...... 57 Occupational safety ...... 12 Open Air (Vinyl) ......24 Operating manuals (types) ...... 5 Operating the lid......43 Optics: ...... 31, 32, 80 P Object holder ...... 51 Occlusion clamp ...... 57

Optional		9
Calibration model		74
Objects (Universal)		69
Occlusion model	55,	59
Single jaw model		46
Tooth stump (individual tooth)		62
Triple Tray® impression		65
Positioning items		69
Positioning objects		69
Power cable		24
Power consumption (technical data)		84
Pre-treatment of dark surfaces		
Pre-treatment of glossy surfaces		45
Protective measures		
R		
Removing dirt		80
Repair		
Representation		
Requirements (workstation)		
Rubber band		
Nabber burid		
S		
Safety	12-	-21
Correct response		
Design		
Injuries		
Material damages		
_		
Symbols		
Touching		
User qualification		
Scanbody		54
Scanning		
Principle		
·		77
Starting with the touchscreen		
Starting with the touchscreen		
Starting with the touchscreen	 35,	81
Starting with the touchscreen	 35,	81
Starting with the touchscreen	 35,	81
Starting with the touchscreen	 35, 35,	81 81
Starting with the touchscreen	35, 35,	81 81
Starting with the touchscreen	35, 35,	81 81 36
Starting with the touchscreen Scope of delivery Sensor number Serial number Setup Dimensions Location	35, 35,	81 81 36 36
Starting with the touchscreen	35, 35, 35,	81 36 36 36
Starting with the touchscreen	35, 35,	81 81 36 36 62 62
Starting with the touchscreen	35, 35,  35,	81 36 36 36 62 62 81
Starting with the touchscreen  Scope of delivery  Sensor number  Serial number  Setup  Dimensions  Location  Weight  Slot (multiDie Adapter)  Slots (multiDie Adapter)  SO number	35, 35,  35, 77,	81 36 36 62 81 79

PC (Scope of delivery)





Squeeze bite	54
Standards (EU Directives)	86
Standards (Safety)	12
Status LED	33
STL (data format)	9
Storage temperature (technical data)	84
Striped pattern	76
Supply voltage (technical data)	84
Symbols	6
System plate	32, 44, 45
Т	
Tick (symbol)	6
Tip (convention)	6
Touchscreen	30, 77
Training	12
Triangle (symbol)	6
Triple Tray®	
Clamp	65
Insert/remove	,
Scope of delivery	
Triple Tray® impression scan	
Type plate	35, 81
U	
Units of measurement	
International System of Units	8

Universal Mode Unpacking	
USB cable	
Use (as intended)	
V	
V (Technical specifications)	84
Vacuuming	80
Vestibular scan (occlusion)	
Vinyl High Resolution	79
W	
Warning (convention)	10
Warranty	81
Watt (technical data)	84
Wax-up	54
WEEE	82
Weight (technical data)	
White light	
Width of the device	
Windows (technical data)	85
Workstation	36





#### Manufacturer

smart optics Sensortechnik GmbH Lise-Meitner-Allee 10 D-44801 Bochum, Germany









ID: SO-20900-OMBA-90334-90351-90400-EN

www.smartoptics.de



