

smart optics scanBox

The scanBox is a fully automated, open 3D dental scanner from smart optics. The scanBox is supplied without CAD software and can be used with every open CAD software.

Alternatively, the scanBox is available in a package with the exocad CAD software - ideal for a start in CAD.

For whom is the scanBox particularly suited?


For every beginner as well as professional users who look to easy, quick and cost-effective manufacturing. The scanBox is a safe investment as there are no licence fees for the scanner.

How about the quality aspect?

All smart optics products are supplied in the usual 'Made in Germany' quality - and this includes the scanBox. A 2-year guarantee rounds things off. The guarantee can be extended to 36 months upon request.

You have questions?

Information on our products is available from our specialized retailers. A comprehensive list of our specialized retailers is available at: www.smartoptics.de.



TECHNICAL DATA

Precision in accordance with DIN EN ISO 12836:2013-01	6 µm
Dimensions (WxHxD)	360 x 310 x 390 mm
Measuring field	80 x 60 x 85 mm
Weight	11 kg
Mains voltage	100 – 240 VAC, 50/60 Hz
Port	USB
Output format	STL



smart optics

smart optics Sensortechnik GmbH
Lise-Meitner-Allee 10 | 44801 Bochum, Germany
Tel.: + 49 234 / 29 828-0 | Fax + 49 234 / 29 828-20
E-mail: info@smartoptics.de | www.smartoptics.de

Ask your specialized retailer for a non-committal demonstration.
smart optics is certified in accordance with DIN EN ISO 13485 and DIN

Made
in
Germany